


## Lab : ✓

# Sample Collection & transport



# General consideration for proper sample collection

- Every laboratory **should provide proper guidelines** for collection of samples.
- All diagnostic information depends on **quality of sample** received.
- If sample collection, transport, media are not proper, it will give **false results**.

- Collect sample **before administering antimicrobial** agents when possible.
- Collect sample with **as little contamination** from indigenous microbiota as possible to ensure that the sample will be representative of the infected site.
- Use **sterile equipment** and **aseptic technique** to collect specimens to prevent introduction of microorganisms during invasive procedures.
- Clearly **label the specimen** container with the patient's name and identification number. Always include date and time of collection.

- Collect an **adequate amount** of specimen. Inadequate amounts of specimen may yield false-negative results.
- Identify the **specimen source** and/or **specific site** correctly so that proper culture media will be selected during processing the laboratory.
- Collect specimens in sturdy, **sterile**, screw-cap, leak proof containers with lids that do not create an aerosol when opened.
- Collect sample after **proper preparation of area** by spirit and povidone iodine.

## Label high risk specimens:

- Sputum with suspected **Tuberculosis**
- Fecal samples suspected with **Cholera, Typhoid.**
- Serum when suspected with **HIV/ HBV/ HCV**, infections



# Rejection criteria

- Leaking/broken container.
- Insufficient amount.
- Improper labeling.
- More time lag between collection and transport.
- Improper transport media.
- Improper transport temperature.
- Hemolyzed sample



## Collection of different samples

- Blood
- Urine
- Stool
- Rectal swab
- CSF

# Collection of blood

- Take proper precaution (gloves).
- Avoid contamination.
- Palpate vein.
- Apply disinfectant.
- Use sterile needle and syringe.
- **Collect appropriate amount of blood depending on the test.**
- Infuse it into tube with anti-coagulants.



# Collection of blood for serology

- Follow the same precautions and method for collections of blood.
- Collect blood into **plain vacutte (test tube)**.
- It will coagulate .
- **Centrifuge** it at 1000-2000 rpm for 5 mins.
- Thus serum will be separated.



# Collection of urine

- Allow first few ml of urine to drain.
- Take mid stream clean catch urine this will prevent contamination.
- Collect it into wide, sterile, leak proof container.
- Transport it directly to the laboratory after proper labeling.
- If delay in transport, preserve it at  $4^{\circ}$  degree centigrade



# Collection of stool

- Collect it into wide mouth, clean, sterile, leak proof container
- Label properly
- Min  $\bullet$  gms is needed if **solid stool** and  $\blacktriangledown$  ml if **liquid** stool.
- It should not be contaminated with urine.
- If not possible to collect it as in children, elderly, debilitated patients, collect rectal swab.
- Do not refrigerate stool.
- If delay in transport preserve it into 10% formalin or buffered glycerol saline.




# Collection of rectal swab

- Take sterile swab.
- Apply it in **anal canal**.
- **Rotate it for 10 seconds**.
- Avoid contact with skin.
- Seal it directly into swab container to prevent contamination.
- **Transport immediately** to lab , if not possible, **preserve at -70 degree** centigrade.




# Collection of cerebrospinal fluid (CSF)

- Collect only 5-10 ml into a labeled sterile container.
- Removal of large volume of CSF lead to headache.
- The fluid to be collected at the rate of 5-10 drops per second.
- Do not refrigerate CSF. If delay is anticipated leave at room temperature.



# Collection of cerebrospinal fluid (CSF)

If sudden removal of fluid is allowed  
may draw down **cerebellum into the  
Foramen magnum and compress the  
Medulla of the Brain.**



# Transportation of samples

- Stool, CSF and sputum should be transported at room temperature, not in refrigerator.
- Urine , swabs, skin samples ,water & food samples are transported a sap at room temp, but if it is not possible, preserve them in refrigerator.